

**AMAZING DISCOUNT ON ANESTHESIA INSTRUMENTS****LARYNGOSCOPES*****IMAGES FOR MECHANICAL EVALUATION*****IMG-0472****MCCOY FLEX TIP SETS**

BLADES: McCoy Flexi Tip blades have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Extremely useful for better visualization of the vocal cords during difficult intubation.

Blades: 0, 1, 2, 3, 4, 5

**IMG-0474****FIBER OPTIC LARYNGOSCOPE SET**

with fiber optic technology that transmits bright, focused light directly to the blade tip. This technology enhances visibility during intubation, making it easier to visualize vocal cords even in low-light conditions. LIBRAN's sets come in various sizes and configurations, including the popular **Macintosh blade design**, compatible with international standards.

Blades: 0, 1, 2, 3, 4, 5

IMG-0473**MCCOY FLEX TIP SETS
INTEGRATED BLADES:**

McCoy Flexi Tip integrated blades have a curved profile similar to the Macintosh but are equipped with a flexible, hinged distal tip for levering. Extremely useful for better visualization of the vocal cords during difficult intubation.

Blades: 1, 2, 3, 4, 5

**IMG-0475****CONVENTIONAL LARYNGOSCOPE SET**

is used for larynx examination, anesthesia, and resuscitation for endotracheal intubation. Macintosh conventional curved blades are available in dull finish hence reduce reflections. This set consist of 6 Macintosh Curved blades sizes 0, 1, 2, 3, 4 & 5 & 1 medium Handle with Conventional Illumination First Responder Kit 100% Brand New Premium Grade Stainless Steel Construction. Built from durable materials that are made to last and serve you for many years.

Blades: 0, 1, 2, 3, 4, 5





IMG-0476

Macintosh Laryngoscope
Fiber Optic – **REUSEABLE**
BLADES: Stainless Steel
HANDLE: Brass C/P.
FINISH: Mirror



IMG-0477

Macintosh Laryngoscope
Fiber Optic (**INTEGRATED-
FLEXI TIP**)– **REUSEABLE**
BLADES: Stainless Steel
HANDLE: Brass C/P.
FINISH: Mirror
Blades: 1, 2, 3, 4, 5



IMG-0478

Flexi-Tip Fiber Optic
Laryngoscope
Blade: 1, 2, 3, 4, 5
Stainless Steel



IMG-0479

FIBER OPTIC
INTEGRATED Flexi-tip
Blades: 1, 2, 3, 4, 5
Stainless Steel



IMG-0480

FIBER OPTIC Flexi-tip
Blades: 1, 2, 3
Stainless Steel



IMG-0481



IMG-0482



IMG-0483



IMG-0484




IMG-0485






IMG-0486

MILLER Laryngoscope Set, CONVENTIONAL,
BLADES #. 1, 2, 3, 4 Stainless Steel.



IMG-0487




IMG-0488

IMG-0489



→LS-010

Macintosh Fiber Optic Laryngoscope Blades
US\$. 5.20 Each



Blade #. 0

Blade #. 1

→LS-012

Miller Port Fiber Optic Blades #. 0 & 1
US\$. 5.30 Each



LIBRAN Surgical

→LS-011

Macintosh Fiber Optic Fiber Optic Blades
#. 00 to 3 US\$. 5.20 Each



→LS-008

FIBER OPTIC BATTERY HANDLE
Large & Med US\$. 5.10
Small, Sty & Mini US\$. 5.00 Each



FLEXI-TIP

→LS-001

Macintosh Fiber Optic Blades
#. 1 to 5 US\$. 5.6 Each



→LS-002

Macintosh Fiber Optic Blades
#. 0 to 3 US\$. 4.90 Each



INTEGRATED

→LS-003

Macintosh Fiber Optic Blades
#. 0 to 3 US\$. 5.50 Each



→LS-004

Macintosh Conventional Blades
#. 00 to 3 US\$. 4.60 Each



→LS-031

Macintosh Conventional Laryngoscope Set of 6 Blades #. 0, 1, 2, 3, 4, 5 Stainless Steel with Medium Battery Handle in Plastic Case.
US\$. 40.25 Each Set.



→LS-022

McCoy Fiber Optic Laryngoscope Set of 3 Pcs. Blade #. 1, 2, 3, Stainless Steel with Medium Battery Handle in Plastic Box.
US\$. 38.60 Each Set.



IMG-0490



**MACINTOSH FIBER OPTIC INTEGRATED FLEXI-TIP LARYNGOSCOPE SET
BLADE #. 1, 2, 3, 4 & 5 STAINLESS STEEL – BATTERY HANDLE MEDIUM BRASS C/P.**

- ▶ Illuminant (xenon lamp) is located in the Battery Handle
- ▶ Cold Light with Flexible Tip For Difficult Intubation.
- ▶ Fiber Optic Handle with 3.5mm² for improved light transmission and longer life.
- ▶ Available in Halogen and LED Illumination Bulbs 2.5V
- ▶ Handles requires 2 x C size batteries (Batteries are not included)

AMAZING DISCOUNT ON ANESTHESIA INSTRUMENTS

LARYNGOSCOPES



→LS-032

Macintosh Conventional Sets, Polished, Stainless Steel, Incl.: Blades Size 0, 1, 2, 3, 4, 5 with Battery Handle Size Medium in Durable Fancy Plastic Box.

US\$. 43.25 Each Set.

~~Reg. Price US\$. 80.10~~



→LS-040

Miller Laryngoscope Fiber Optic Set of 5 Blades, Stainless Steel, incl.: Blades #.0, 1, 2, 3, 4, 5 with Medium Battery Handle in Durable fancy Plastic Box.

US\$. 47.50 Each Set.

~~Reg. Price US\$. 75.10~~



→LS-039

Miller Laryngoscope Fiber Optic Set of 4 Blades Stainless Steel, Blade #.1, 2, 3, 4 with Medium Battery Handle in Pouch/Bag.

US\$. 33.10

~~Reg. Price US\$. 62.30~~



→LS-033

Macintosh Laryngoscope Fiber, Optic, Flexi-tip, Stainless Steel, Blades Size 1, 2, 3, 4, 5 with Medium Battery Handle, in Fancy Plastic Case.

US\$. 45.10 Each Set.

~~Reg. Price US\$. 65.25~~



→LS-034

Macintosh Laryngoscope Fiber Optic Flexi-Tip, Integrated Sets, Stainless Steel Blade #. 1, 2, 3, 4, 5 Battery Handle Medium, in Fancy Plastic Box

US\$. 48.55 Each Set.

~~Reg. Price US\$. 67.00~~



→LS-041

Macintosh Laryngoscope Set of 5 Blades #.1,2,3,4,5 Stainless Steel, Fiber Optic, Battery Handle Medium in Plastic Case.


US\$. 56.90 Each Set.

~~Reg. Price US\$. 62.30~~

AMAZING DISCOUNT ON ANESTHESIA INSTRUMENTS

LARYNGOSCOPES

→LS-050



→LS-051



Macintosh Mega (Integrated) Fiber Optic Laryngoscope Set of 5 Blade #.1, 2, 3, 4 & 5 Stainless Steel Blades + Any ONE Battery Handle Large, Medium, Penlite Stubby in Plastic Case.

US\$. 65.10 Each Sets.

REG. PRICE: US\$. 72.10

Macintosh Laryngoscope American Profile Set, Set of 6 Conventional, Blade No. 0, 1, 2, 3, 4, 5 Stainless Steel, Battery Handle Medium Brass Chrome Plated (OR) in Stainless Steel without extra charges, Set in Fancy durable Plastic Case.

US\$. 77.00 Each Set.

REG. PRICE: US\$. 82.50

IMAGES FOR MECHANICAL EVALUATION

**Medical LED Light Otoscope Set For Veterinary,
ENT Fiber Optic Magnifying Lens, Brass
Chrome Plated.**



REG. PRICE US\$ 68.50 PRICE: US\$. 52.00



REG. PRICE LS-051 to LS-053 US\$ 79.20 Each

PRICE LS-051 to LS-053 US\$. 68.10 Each



PRICE US\$. 49.00



PRICE US\$. 49.00



I S-056

PRICE US\$.42.00



I S-057

PRICE US\$.52.00



I S-058

PRICE US\$.58.20

AMAZING DISCOUNT ON ANESTHESIA INSTRUMENTS

OTOSCOPES AND OPHTHALMOSCOPES



OTOSCOPE+OPHTHALMOSCOPE
FIBER OPTIC
INCLUDING
 Ear Funnels: Ø2.5/Ø4.0mm -10 x 10
 In Fancy Rexene Pouch.
 US\$. 33.60 Each Set

OTOSCOPE - FIBER OPTIC
INCLUDING EAR FUNNELS:
 8 Pcs. of Ø2.4mm for Infants - 8 Pcs. of 3mm for Child
 8 Pcs. Of 4mm for Adult - 8 Pcs. Of 5mm for Adult Large
 In Fancy Rexene Pouch.
 US\$. 33.60 Each Set



1x OTOSCOPE + 1x OPHTHALMOSCOPE SET (Blue), INCLUDES 10 Specula of 2mm & 10 specula 4mm in 1 Silver Carrying Case. –
 Perfect for Nursing & Medical Students with Carry Case. (AVAILABLE IN ANY COLOR)
 US\$. 60.60 Each Set.

AMAZING DISCOUNT ON ANESTHESIA INSTRUMENTS

OTOSCOPES AND OPHTHALMOSCOPES

→LS-123



DIFFERENT COLORS

Mini Fiber Optic Otoscope (Pocket Size) Medical Ent Diagnostic Set, with 5X 2.5mm/5x 4mm disposable plastic specula
US\$. 16.50 Each Set.

→LS-124



Mini Fiber Optic Otoscope (Pocket Size) Medical Ent Diagnostic Set, with 5X 2.5mm/5x 4mm disposable plastic specula

US\$. 16.50 Each Set.

→LS-125



Premium Fiber Optic Otoscope Ophthalmoscope Examination LED Diagnostic ENT SET. (BLACK)
10X Disposable Plastic Specula in durable fancy Plastic case.

US\$. 56.10 Each Sets.

→LS-126



Premium Fiber Optic Otoscope Ophthalmoscope Examination LED Diagnostic ENT SET. (BLUE)
10X Disposable Plastic Specula in durable fancy Plastic case.

US\$. 46.10 Each Sets.

→LS-127



Otoscope/Ophthalmoscope Set with 9+ Reusable Cannula
In durable fancy Box.

US\$. 39.90 Each Set.

→LS-128



REUSEABLE

LIBRAN Quality Ophthalmoscope For Professionals and Students Quality Tool Examine Eyes. Packed in Top Quality Fancy/Durable Plastic Case.

US\$. 22.10 Each Set.



LS-129

**FIBER OPTIC Oscope +
Ophthalmoscope Sets in Fancy
Durable Box.**

~~Regular Price US\$ 55.50~~

PRICE US\$ 39.25 Each Set.



**FIBER OPTIC
OTOSCOPE**
10 plastic cannula

LS-130



**FIBER OPTIC
OTOSCOPE**
10 Plastic Cannula

LS-131



Wide-angle viewing lens
allows for instrumentation
under magnification

LS-132



Mini Pocket Oscope
Set LED Ear Scope with 5
Tips (Black)

LS-133



LS-134

**Otoscope Set Durable Build
with Bright LED for ENT
Surgical Use, A+ Grade Quality,**

PRICE US\$ 40.20



LS-135

**Absolutely Great and Most Advanced
Quality Fiber Optic Oscope. Truly
provides for the best diagnosis.**

PRICE US\$ 52.20



PRICE US\$ 43.20

FIBER-OPTIC ILLUMINATION

POWER BUTTON ON/OFF

SPECULA JACK

3 X MAGNIFIER

IMAGES FOR MECHANICAL EVALUATION



IMG-0620



IMG-0621



IMG-0622



IMG-0623



IMG-0624



IMG-0625



IMG-0626



IMG-0627



IMG-0628

REUSABLE EAR SPECULA STAINLESS STEEL



IMG-0629

AMAZING DISCOUNT ON GYNECOLOGY INSTRUMENTS



SPECULUMS

CUSCO



→LS-220

Size	Beak Length	Width
Ex-Small	ca. 7.5cm	ca. 1.7cm
Small	ca. 8.0cm	ca. 2.2cm
Medium	ca. 10.0cm	ca. 2.5cm
Large	ca. 10.5cm	ca. 2.7cm

Made from Stainless Steel AISI-420, used for visual inspection of the Vagina.
Available in Satin & Bright Mirror Finish.

US\$. 7.45 Each Piece.

COLLIN



→LS-221

Small: 85 x 30 mm to 90 x 29 mm
Medium: 100 x 35 mm
Large: 110 x 40 mm to 110 x 43 mm
Extra Small/Specialty: 65 x 20 mm
Narrow/Specific: 90 x 16 mm
Fig. 1 (Small): 90 x 23-25 mm
Made from Stainless Steel AISI-420, used for visual inspection of the Vagina.
Available in Satin & Bright Mirror Finish.
US\$. 10.95 Each Piece.

PEDERSON



→LS-222

SMALL: 7.5cm x 1.5cm
MMEDIUM: 10.0cm x 2.5cm
LARGE: 10.0cm x 2.5cm
X-LARGE: 11.5cm x 2.5cm
X-NARROW: 11.5cm x 1.3cm
XXL: 18.0cm x 2.5cm
Made from Stainless Steel AISI-420, used for visual inspection of the Vagina.
Available in Satin & Bright Mirror Finish.
US\$. 35.95 Each Piece.

GRAVES



→LS-223

Extra Small: 3" x 3/4"
Small: 3-1/2" Xx7/8"
Medium Size: 1-3/8" x 4" Angled
Medium: 4"x1- 1/4"
Large: 4-1/2" x 1 3/8"
Extra Large: 5" x 1 1/2"
Made from Stainless Steel AISI-420, used for visual inspection of the Vagina.
Available in Satin & Bright Mirror Finish.
US\$. 18.95 Each Piece.

For CATTLES



→LS-224

Working Length: 33cm
Maximum Open: 11cm
-Round head design, protect inner wall of cervical
-Easy clean, high temperature disinfection can be performed
Made from Stainless Steel AISI-420, used for visual inspection of the Vagina.
Available in Satin & Bright Mirror Finish.
US\$. 38.90 Each Piece.

LB NEW SPEC.



→LS-225

New designed Vaginal Speculums, Small, Medium and Large.
Made from Stainless Steel AISI-420, used for visual inspection of the Vagina.
Available in Satin & Bright Mirror Finish.
US\$. 38.90 Each Piece.

AMAZING DISCOUNT ON GYNECOLOGY INSTRUMENTS

SPECULUMS

MILLER



→LS-226


Miller Vaginal Speculum,
Size: Working Length 10.1cm

Made from Stainless Steel AISI-420, used
for visual inspection of the Vagina.

Available in Satin & Bright Mirror Finish.

US\$. 22.95 Each Piece.

GUTTMANN



→LS-227

Guttman Vaginal Speculum
have a range of combinations
such as speculums with self-
retaining, Examining type, Size
1-3/8" x 4".

Made from Stainless Steel AISI-420, used
for visual inspection of the Vagina.

Available in Satin & Bright Mirror Finish.

US\$. 22.95 Each Piece.

WINTERTON (CUSCO)



→LS-228

WINTERTON VAGINAL SPECULUM
LARGE. BLADE 25MM x 115MM.

Made from Stainless Steel AISI-420, used
for visual inspection of the Vagina.

Available in Satin & Bright Mirror Finish.

US\$. 29.95 Each Piece.




→LS-229

SIMS VAGINAL SPECULUM
Small, Medium, Large, Extra
Large.

Made from Stainless Steel
Available in Satin & Polish

US\$. 3.50 Each Piece.



→LS-230

SIMS VAGINAL SPECULUMS

65 x 25mm, 70 x 30mm (small)

70 x 30mm, 70 x 35mm (medium)

85 x 35mm, 90 x 40mm (large)

•Smooth, Curved Blade

• Long handle for a comfortable grip

• Lightweight design

• Sterilizable and reusable

Made from Stainless Steel

Available in Satin & Polish

US\$. 3.50 Each Piece.



32 x 30mm

60 x 25mm

→LS-232

SIMS vaginal speculum

25 x 60 mm, 30 x 60 mm, 12.5 cm.

Made from Stainless Steel

Available in Satin & Polish

US\$. 3.50 Each Piece.



→LS-231

SIMS VAGINAL SPECULUMS

SINGLE END

SIZES

PRICE

▶ 70mm x 27mm **US\$. 6.70 Each.**

▶ 80mm x 30mm **US\$. 6.80 Each.**

▶ 90mm x 36mm **US\$. 6.90 Each.**

▶ 110mm x 27mm **US\$. 7.00 Each.**

▶ 110mm x 36mm **US\$. 7.10 Each.**

Made from Stainless Steel AISI-410

Available in Satin & Polish

AMAZING DISCOUNT ON GYNECOLOGY INSTRUMENTS



AUVARD WEIGHTED VAGINAL SPECULUMS

→LS-232



AUVARD WEIGHTED VAGINAL SPECULA

Slightly Curved Blade 5.5" x 1.25", at 90° Angle.
Weight: 0.59 LBS.
Single Ended.
Made from Stainless Steel.
US\$. 19.20 Each Piece.

→LS-233



AUVARD WEIGHTED VAGINAL SPECULA

2.0lbs. - 7.0cm x 4.0cm
US\$. 23.50 Each Piece
2.5lbs. - 7.0cm x 4.0cm
US\$. 25.50 Each Piece
2.5lbs. - 10.0cm x 4.5cm
US\$. 28.50 Each Piece
3.0lbs. - 7.0cm x 4.0cm
US\$. 27.50 Each Piece
Made from Stainless Steel.

→LS-234



AUVARD WEIGHTED VAGINAL SPECULA

80 x 43cm
Length: 23.5cm
Weight: 2.5Lbs.
Made from Stainless Steel.
US\$. 23.20 Each Piece

→LS-235



AUVARD WEIGHTED VAGINAL SPECULA

230 mm (9"),
Working Length: 65 mm,
Jaw Width 28 mm-38 mm,
Weight: 1,530 g
Made from Stainless Steel.
US\$. 27.20 Each Piece



AUVARD WEIGHTED VAGINAL SPECULA

230 mm (9"), work. length: 65 mm, jaw width 28 mm - 38 mm, weight: 1,530 g
Made from Stainless Steel.
US\$. 13.20 Each Piece.



AUVARD WEIGHTED VAGINAL SPECULA

a-2.0lbs, 75°-Blade 38mm x 83mm
b-2.5lbs, 75°-Blade 38mm x 83mm
c-3.0lbs, 75°-Blade 38mm x 83mm
d-2.5lbs, 45°-Blade 35mm x 102mm
e-2.5lbs, 90°-Blade 25mm x 140mm

Made from Stainless Steel.

a-US\$. 13.15 Each Piece.

b-US\$. 14.15 Each Piece.

c-US\$. 15.30 Each Piece.

d-US\$. 14.40 Each Piece.

e-US\$. 15.30 Each Piece.



AUVARD WEIGHTED VAGINAL SPECULA

Used in gynecological procedures and it features a blade that is bent at a right angle with a spherical weight on the handle.

Broad Jaw for Better Visualization.

Made from German grade Surgical Stainless Steel.
Matte Finish,
Non-sterile and reusable

US\$. 16.15 Each Piece.



AUVARD WEIGHTED VAGINAL SPECULA

Overall Length: 9 inch
Weight: 3 Lbs.

Made from Stainless Steel.

US\$. 13.20 Each Piece.



AUVARD WEIGHTED VAGINAL SPECULA

140mm X 30mm Blade
Goose Neck
Weighted 2.5 lbs
Made from Stainless Steel.
US\$. 12.15 Each Piece.

